

OTICON | More

life-changing  
technology

# More potential than ever

In a wide range of innovative styles

## NEW Innovative miniBTE styles

Choosing Oticon More™ has never been easier. With two new miniBTE styles joining our flagship miniRITE styles, you can now bring the unique sound quality of Oticon More to more clients who need or want a BTE hearing aid.

An entire family of hearing aids now stands ready to face diverse client needs and preferences. And ready to give more people the access to sound they need.



## More to offer altogether:

- **NEW miniBTE R** style suitable for a hook and mould or a thin tube and dome, with rechargeable batteries and a SmartCharger or desk charger
- **NEW miniBTE T** style suitable for a hook and mould or a thin tube and dome, with disposable size 312 batteries

**The flexibility of the new miniBTE styles** helps you make the perfect choice for each client with the possibility of using either a hook and mould or a thin tube and dome, and no wax filters to change.

- **Discreet miniRITE R** style featuring a receiver in the ear, with rechargeable batteries and a SmartCharger or desk charger
- **Discreet miniRITE T** style featuring a receiver in the ear and disposable size 312 batteries


## All rechargeable Oticon More styles

give your clients an appealing choice of power sources. Let them choose between a portable and sleek SmartCharger, or a desk charger.




# Covering a wide range of hearing loss

## miniRITE



## miniBTE



# Outstanding features at three price points

Benefits	Features	Oticon More 1	Oticon More 2	Oticon More 3
Speech Understanding	MoreSound Intelligence	Level 1	Level 2	Level 3
	- Environment configuration	5 Options	5 Options	3 Options
	- Virtual Outer Ear	3 Configurations	1 Configuration	1 Configuration
	- Spatial Balancer	100%	60%	60%
	- Neural Noise Suppression, Difficult/Easy	10 dB/4 dB	6 dB/2 dB	6 dB/0 dB
	- Sound Enhancer	3 Configurations	2 Configurations	1 Configuration
	MoreSound Amplifier	•	•	•
	Feedback Prevention	MoreSound Optimizer & Feedback shield	MoreSound Optimizer & Feedback shield	MoreSound Optimizer & Feedback shield
	Spatial Sound	4 Estimators	2 Estimators	2 Estimators
	Soft Speech Booster	•	•	•
Frequency lowering	Speech Rescue	Speech Rescue	Speech Rescue	
Sound Quality	Clear Dynamics	•	•	-
	Better-Ear Priority	•	•	-
	Fitting Bandwidth*	10 kHz	8 kHz	8 kHz
Listening Comfort	Transient Noise Management	4 configurations	3 configurations	3 configurations
	Wind Noise Management	•	•	•
Personalization and Optimizing Fitting	Fitting Bands	24	20	18
	Multiple Directionality options	•	•	•
	Adaptation Manager	•	•	•
	Fitting Formulas	VAC+, NAL-NL1/NAL-NL2, DSL 5.0	VAC+, NAL-NL1/NAL-NL2, DSL 5.0	VAC+, NAL-NL1/NAL-NL2, DSL 5.0
	Tinnitus SoundSupport	•	•	•
	CROS/BiCROS support	•	•	•
Connecting to the world	Direct streaming**	•	•	•
	Oticon ON and Oticon RemoteCare apps	•	•	•
	ConnectClip and EduMic	•	•	•
	TV Adapter 3.0 and Phone Adapter 2.0	•	•	•
	Remote Control 3.0	•	•	•

\* Bandwidth accessible for gain adjustments during fitting \*\* From iPhone®, iPad®, iPod touch®, and selected Android™ devices